

MX-S+


Digital Vortex Mixer

Features

- LCD screen displays real-time speed and time for easy adjustment and monitoring.
- Adjustable speed ranges from gentle to vigorous mixing to suit various needs.
- Stable with pressure-mounted accessories to prevent slippage and splashing.
- Brushless DC motor for high stability and maintenance-free operation.
- Sturdy base ensures stability during operation for safety.
- Includes various accessories for different sizes and types of tubes and plates.
- Low noise for a quiet working environment.



Specifications	MX-S+
Speed	100-3000 rpm
Display	LCD
Amplitude	4 mm
Operating Modes	Touch, Continuous
Time Range	1 s - 99 min 59 s
Motor Type	Brushless DC Motor
Dimensions (L×W×H)	160x146x115mm
Weight	3.4kg
Power Supply	100-240 V, 50/60 Hz
Input Voltage	24 V
Power Consumption	25 W
Ambient Temperature & Humidity	5-40°C, 80% RH
Enclosure Protection	IP21



Accessories of Vortex Mixer MX-S+ & MX-Pro

Accessories of MX-S+&MX-Pro

VT2.3 Cat. No.18901365

Universal top plate, compatible with VT1.3.1-VT1.3.12 Tube Adapter



VT2.5 Cat. No.18901367

Erlenmeyer flask attachment (100mL-250mL), Operating speed: 0~1500rpm



VT2.4 Cat. No.18901366

Ampule/Tube Attachment, 4pc ampule/tube (0.2mL-15mL)



PS2.1 Cat. No.18901360

Microplate clamp, Suitable for workboards (full skirt PCR plates, cell culture plates, enzyme-linked immunosorbent assay plates, deep well plates, etc. Half skirt/no skirt PCR plates need to be used with rectangular foam, which is already included)



Accessories of MX-S+&MX-Pro (New)

VT2.1 Cat. No.18901379

Top cup, Suitable for MX-Pro



VT2.6.2 Cat. No.18901377

Tube Adapter, 21 holes, Ø8mm(1.5-2mL), Applicable to RUN mode



VT2.1 Cat. No.18901380

Top cup, Suitable for MX-S+



VT2.6.3 Cat. No.18901376

Tube Adapter, 16 holes, Ø10mm(2-5mL, Not compatible with 5 mL conical tubes. Compatible with standard 5 mL centrifuge tubes.), Applicable to RUN mode



VT2.2 Cat. No.18901373

Flat Rubber Sleeve (Standard)



VT2.6.4 Cat. No.18901375

Tube Adapter, 8 holes, Ø14mm (10-15mL, supports 5 mL conical tubes; incompatible with standard 5 mL centrifuge tubes.), Applicable to RUN mode



VT2.6 Cat. No.18901368

Snap-in Universal top plate, compatible with VT2.6.1-VT2.6.5 Tube Adapter



VT2.6.5 Cat. No.18901374

Tube Adapter, 4 holes, Ø25.5mm(50mL), Applicable to RUN mode



VT2.6.1 Cat. No.18901378

Tube Adapter, 25 holes, Ø4.5mm(0.1-0.5mL), Applicable to RUN mode



*Material color may vary slightly during the initial launch phase. Please refer to the actual product.

Accessories of MX-S+&MX-Pro (New)

VT2.6.6	Cat. No.18902244	Tube Adapter, 19 holes × Ø9mm, Applicable to 2 mL cryovial, Applicable to RUN mode 	VT 2.7.5	Cat. No.18901871	Horizontal Plastic Clamp, 12 holes × 15mL Test Tube (Circular plate) 
VT2.6.7	Cat. No.18902280	Vial adapter, 8 holes, Ø17mm (suitable for 5 mL vials), Applicable to RUN mode 	VT 2.7.6	Cat. No.18901872	Horizontal Plastic Clamp, 6 holes × 50mL Test Tube (Circular plate) 
VT2.6.8	Cat. No.18902281	Vial adapter, 6 holes, Ø20.5mm (suitable for 10 mL vials), Applicable to RUN mode 	VT 2.8.1	Cat. No.18902336	50 mL vertical multi-format rack, 12 x 50 mL, including 24 clamps, suitable for RUN mode 
VT2.6.9	Cat. No.18902282	Vial adapter, 4 holes, Ø26.5mm (suitable for 20 mL vials), Applicable to RUN mode 	VT 2.8.2	Cat. No.18902337	15 mL vertical multi-format rack, 12 x 15mL, including 24 clamps, suitable for RUN mode 
VT 2.6.10	Cat. No.18902283	Vial adapter, 4 holes, Ø26.5mm (suitable for 30 mL vials), Applicable to RUN mode 	VT 2.8.3	Cat. No.18902338	1.5 mL vertical multi-format rack, 12 x 1.5mL, including 24 clamps, suitable for RUN mode 
VT 2.7.4	Cat. No.18901870	Horizontal Plastic Clamp, 24 holes × 1.5mL Test Tube (Circular plate) 	PS2.2	Cat. No.18901361	Snap-in Microplate clamp, Suitable for workboards (full skirt PCR plates, cell culture plates, enzyme-linked immunosorbent assay plates, deep well plates, etc. Half skirt/no skirt PCR plates need to be used with rectangular foam, which is already included) 

*Material color may vary slightly during the initial launch phase. Please refer to the actual product.